

B3GALT4 Antibody



PACO22195

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Involved in GM1/GD1B/GA1 ganglioside biosynthesis.

Gene ID:

B3GALT4

Uniprot

O96024

Synonyms:

b3Gal-T4; B3GT4; beta-1; 3-galactosyltransferase 4; beta-1


Immunogen:

Synthesized peptide derived from internal of human B3GALT4.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from Jurkat cells, using B3GALT4 antibody.