

BMP3 Antibody



PACO22205

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Negatively regulates bone density. Antagonizes the ability of certain osteogenic BMPs to induce osteoprogenitor differentiation and ossification.

Gene ID:

BMP3

Uniprot

P12645

Synonyms:

Bone morphogenetic protein 3; BMP-3; Osteogenin; BMP-3A; Bone morphogenetic protein 3A


Immunogen:

Synthesized peptide derived from internal of human BMP3.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HeLa cells, HuvEc cells and HepG2 cells, using BMP3 antibody.