

CA3 Antibody



PACO22233

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Reversible hydration of carbon dioxide.

Gene ID:

CA3

Uniprot

P07451

Synonyms:

CA-III; CAH3; carbonate dehydratase III; carbonic anhydrase 3; carbonic anhydrase III


Immunogen:

Synthesized peptide derived from internal of human CA3.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HUVEC cells, using CA3 antibody.