
Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Delivers copper to copper zinc superoxide dismutase (SOD1).

Gene ID:

CCS

Uniprot

O14618

Synonyms:

copper chaperone for superoxide dismutase; MGC138260;


Immunogen:

Synthesized peptide derived from C-terminal of human CCS.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from Jurkat cells and 3T3 cells, using CCS antibody.