

PLA2G4C Antibody



PACO22300

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Has a preference for arachidonic acid at the sn-2 position of phosphatidylcholine as compared with palmitic acid.

Gene ID:

PLA2G4C

Uniprot

Q9UP65

Synonyms:

cPLA2-gamma; cytosolic phospholipase A2 gamma precursor; phospholipase A2 group IVC;


Immunogen:

Synthesized peptide derived from internal of human PLA2G4C.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from Jurkat cells, using PLA2G4C antibody.