

GPRIN1 Antibody



PACO22339

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

May be involved in neurite outgrowth By similarity.

Gene ID:

GPRIN1

Uniprot

Q7Z2K8

Synonyms:

G protein-regulated inducer of neurite outgrowth 1; KIAA1893;


Immunogen:

Synthesized peptide derived from internal of human GPRIN1.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HT-29 cells, using GPRIN1 antibody.