

HLA-H Antibody



PACO22369

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Involved in the presentation of foreign antigens to the immune system.

Gene ID:

HLA-H

Uniprot

P01893

Synonyms:

HLA class I histocompatibility antigen; alpha chain H precursor; HLA-12.4; HLA-AR; HLAH


Immunogen:

Synthesized peptide derived from internal of human HLAH.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from LOVO cells, using HLAH antibody.