

NKX2-4 Antibody



PACO22378

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Probable transcription factor.

Gene ID:

NKX2-4

Uniprot

Q9H2Z4

Synonyms:

Homeobox protein Nkx-2.4; Homeobox protein NKX2.4; Homeobox protein NK-2 homolog D;


Immunogen:

Synthesized peptide derived from internal of human NKX24.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HeLa and COLO205 cells, using NKX24 antibody.