

# KLHL29 Antibody



PACO22382

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:50-1:100

**Protein Background:**

**Gene ID:**

KLHL29

**Uniprot**

Q96CT2

**Synonyms:**

Kelch-like protein 29; Kelch repeat and BTB domain-containing protein 9;

**Immunogen:**

Synthesized peptide derived from internal of human KLHL29.

**Storage:**

Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---

Immunohistochemistry analysis of paraffin-embedded human brain tissue using KLHL29 antibody.

