

PACO22384

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Catalyzes the post-translational modification of target proteins by phosphopantetheine. Can transfer the 4'-phosphopantetheine moiety from coenzyme A to a serine residue of a broad range of acceptors, such as the acyl carrier domain of FASN.

Gene ID:

AASDHPPT

Uniprot

Q9NRN7

Synonyms:

4-phosphopantetheinyl transferase; AASD-PPT; ADPPT; alpha-aminoadipic semialdehyde dehydrogenase-phosphopantetheinyl transferase; EC 2.7.8. -


Immunogen:

Synthesized peptide derived from internal of human AASDHPPT.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from COLO cells, using AASDHPPT antibody.