

PACO22394

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:100-1:500

Protein Background:

May influence APP processing, resulting in a decrease in sAPP-alpha production and increased amyloidogenic P3 peptide production By similarity.

Gene ID:

LDLRAD3

Uniprot

Q86YD5

Synonyms:

Low-density lipoprotein receptor class A domain-containing protein 3;

Immunogen:

Synthesized peptide derived from internal of human LDLRAD3.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

Immunofluorescence analysis of A549 cells, using LDLRAD3 antibody.

