

PACO22408

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:50-1:100

Protein Background:

Elicits growth inhibition on melanoma cells in vitro as well as some other neuroectodermal tumors, including gliomas.

Gene ID:

MIA

Uniprot

Q16674

Synonyms:

Melanoma-derived growth regulatory protein; Melanoma inhibitory activity;

Immunogen:

Synthesized peptide derived from C-terminal of human MIA.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

Immunohistochemistry analysis of paraffin-embedded human heart tissue using MIA antibody.

