

KIF20B Antibody



PACO22429

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Plus-end-directed motor enzyme that is required for completion of cytokinesis.

Gene ID:

KIF20B

Uniprot

Q96Q89

Synonyms:

M-phase phosphoprotein 9; MPP9;


Immunogen:

Synthesized peptide derived from C-terminal of human MPHOSPH9.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from A549 cells, using MPHOSPH9 antibody.