

NPHP4 Antibody



PACO22456

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:50-1:100

Protein Background:

Involved in the organization of apical junctions in kidney cells together with NPHP1 and RPGRIP1L/NPHP8 By similarity. Does not seem to be strictly required for ciliogenesis By similarity.

Gene ID:

NPHP4

Uniprot

O75161

Synonyms:

Nephrocystin-4; Nephroretinin;

Immunogen:

Synthesized peptide derived from internal of human NPHP4.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

Immunohistochemistry analysis of paraffin-embedded human colon carcinoma tissue using NPHP4 antibody.

