

URB1 Antibody



PACO22484

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:50-1:100

Protein Background:

Gene ID:

URB1

Uniprot

O60287

Synonyms:

Nucleolar pre-ribosomal-associated protein 1; Protein URB1;

Immunogen:

Synthesized peptide derived from C-terminal of human URB1.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

Immunohistochemistry analysis of paraffin-embedded human brain tissue using URB1 antibody.

