

PDXDC1 Antibody



PACO22533

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Belongs to the group II decarboxylase family. 2 isoforms of the human protein are produced by alternative splicing.

Gene ID:

PDXDC1

Uniprot

Q6P996

Synonyms:

EC 4.1.1. -; KIAA0251; LP8165; PDXD1; pyridoxal-dependent decarboxylase domain containing 1


Immunogen:

Synthesized peptide derived from internal of human PDXDC1.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HepG2 cells and HUVEC cells, using PDXDC1 antibody.