

PDK2 Antibody



PACO22535

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Gene ID:

PDK2

Uniprot

Q15119

Synonyms:

EC 2.7.11.2; kinase Pyruvate dehydrogenase kinase 2; PDK2; Pyruvate dehydrogenase [lipoamide] kinase isozyme 2; mitochondrial precursor


Immunogen:

Synthesized peptide derived from internal of human PDK2.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HepG2 cells, using PDK2 antibody.