

RAB2B Antibody



PACO22553

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Required for protein transport from the endoplasmic reticulum to the Golgi complex Potential.

Gene ID:

RAB2B

Uniprot

Q8WUD1

Synonyms:

Ras-related protein Rab-2B;


Immunogen:

Synthesized peptide derived from internal of human RAB2B.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from K562 cells, HeLa cells and 293 cells, using RAB2B antibody.