

# RAB29 Antibody



PACO22558

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:3000

**Protein Background:**

**Gene ID:**

RAB29

**Uniprot**

O14966

**Synonyms:**

RAB7; member RAS oncogene family-like 1; RAB7L;

**Immunogen:**


Synthesized peptide derived from internal of human RAB7L1.

**Storage:**

Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---



Western blot analysis of extracts from HT-29 cells, using RAB7L1 antibody.