

ARHGEF12 Antibody



PACO22572

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

May play a role in the regulation of RhoA GTPase by guanine nucleotide-binding alpha-12 (GNA12) and alpha-13 (GNA13). Acts as guanine nucleotide exchange factor (GEF) for RhoA GTPase and may act as GTPase-activating protein (GAP) for GNA12 and GNA13.

Gene ID:

ARHGEF12

Uniprot

Q9NZN5

Synonyms:

ARHGC; ARHGEF12; Leukemia-associated RhoGEF; RHO guanine nucleotide exchange factor (GEF) 12 (LARG); Rho guanine nucleotide exchange factor 12


Immunogen:

Synthesized peptide derived from internal of human ARHGEF12.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from RAW264.7 cells, using ARHGGEF12 antibody.