

ARHGEF19 Antibody



PACO22574

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:50-1:100

Protein Background:

Acts as guanine nucleotide exchange factor (GEF) for RhoA GTPase By similarity.

Gene ID:

ARHGEF19

Uniprot

Q8IW93

Synonyms:

Rho guanine nucleotide exchange factor 19; ARHGEF19;

Immunogen:

Synthesized peptide derived from internal of human ARHGEF19.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue using ARHGEF19 antibody.

