

RRS1 Antibody



PACO22580

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Involved in ribosome biogenesis By similarity.

Gene ID:

RRS1

Uniprot

Q15050

Synonyms:

ribosome biogenesis regulatory protein; RRR;


Immunogen:

Synthesized peptide derived from internal of human RRS1.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from 293 cells and HuvEc cells, using RRS1 antibody.