

REXO1 Antibody



PACO22582

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:50-1:100

Protein Background:

Seems to have no detectable effect on transcription elongation in vitro.

Gene ID:

REXO1

Uniprot

Q8N1G1

Synonyms:

EC 3.1. -. -; Elongin A-binding protein 1; EloA-BP1; Transcription elongation factor B polypeptide 3-binding protein 1;

Immunogen:

Synthesized peptide derived from internal of human REXO1.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

Immunohistochemistry analysis of paraffin-embedded human brain tissue using REXO1 antibody.

