

SH3GL2 Antibody



PACO22612

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Implicated in synaptic vesicle endocytosis. May recruit other proteins to membranes with high curvature.

Gene ID:

SH3GL2

Uniprot

Q99962

Synonyms:

Endophilin-1; SH3 domain-containing GRB2-like protein 2; SH3 domain protein 2A; EEN-B1;


Immunogen:

Synthesized peptide derived from internal of human SH3GL2.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from K562 cells and HeLa cells, using SH3GL2 antibody.