

ATP1B3 Antibody



PACO22622

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

This is the non-catalytic component of the active enzyme, which catalyzes the hydrolysis of ATP coupled with the exchange of Na⁺ and K⁺ ions across the plasma membrane. The exact function of the beta-3 subunit is not known.

Gene ID:

ATP1B3

Uniprot

P54709

Synonyms:

AT1B3; ATP1B3; ATPase; Na⁺/K⁺ transporting; beta 3 polypeptide


Immunogen:

Synthesized peptide derived from Internal of human ATP1B3.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from K562 cells, using ATP1B3 antibody.