

TBX18 Antibody



PACO22644

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Probable transcriptional regulator involved in developmental processes.

Gene ID:

TBX18

Uniprot

O95935

Synonyms:

T-box protein 18; T-box transcription factor TBX18; TBX18;


Immunogen:

Synthesized peptide derived from internal of human TBX18.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from Jurkat cells, using TBX18 antibody.