

USP30 Antibody



PACO22668

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

May participate in the maintenance of mitochondrial morphology.

Gene ID:

USP30

Uniprot

Q70CQ3

Synonyms:

flj40511; mgc10702; ubiquitin specific peptidase 30; ubiquitin-specific proteinase 30; ubp30


Immunogen:

Synthesized peptide derived from internal of human USP30.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HeLa cells, using USP30 antibody.