

USP45 Antibody



PACO22673

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Belongs to the peptidase C19 family. 3 isoforms of the human protein are produced by alternative splicing.

Gene ID:

USP45

Uniprot

Q70EL2

Synonyms:

Ubiquitin carboxyl-terminal hydrolase 45; Ubiquitin thioesterase 45; Ubiquitin-specific-processing protease 45; Deubiquitinating enzyme 45; USP45


Immunogen:

Synthesized peptide derived from internal of human USP45.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from A549 cells, using USP45 antibody.