

UBFD1 Antibody



PACO22675

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

May play a role as NF-kappa-B regulator.

Gene ID:

UBFD1

Uniprot

O14562

Synonyms:

FLJ38870; FLJ42145; ubiquitin family domain containing 1; UBPH;


Immunogen:

Synthesized peptide derived from internal of human UBFD1.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from K562 cells, using UBFD1 antibody.