

UBAC1 Antibody



PACO22677

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Non-catalytic subunit of the KPC complex that acts as E3 ubiquitin-protein ligase. Required for poly-ubiquitination and proteasome-mediated degradation of CDKN1B during G1 phase of the cell cycle.

Gene ID:

UBAC1

Uniprot

Q9BSL1

Synonyms:

E3 ubiquitin-protein ligase KPC2; glioblastoma cell differentiation-related protein 1; Kip1 ubiquitination-promoting complex 2; UBAC1; UBADC1


Immunogen:

Synthesized peptide derived from internal of human UBAC1.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from 293 cells and HUVEC cells, using UBAC1 antibody.