

ZNF771 Antibody



PACO22696

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

May be involved in transcriptional regulation.

Gene ID:

ZNF771

Uniprot

Q7L3S4

Synonyms:

Zinc finger protein 771; Mesenchymal stem cell protein DSC43; ZNF771;


Immunogen:

Synthesized peptide derived from internal of human ZNF771.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HepG2 cells, using ZNF771 antibody.