

RNF113B Antibody



PACO22721

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:**Gene ID:**

RNF113B

Uniprot

Q8IZP6

Synonyms:

bA10G5.1; MGC26599; R113B; RING finger protein 113B; RNF113B


Immunogen:

Synthesized peptide derived from Internal of human RNF113B.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from COLO205/K562/Jurkat/293 cells, using RNF113B antibody.