

# CYP11B1/CYP11B2 Antibody



PACO22740

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:3000

**Protein Background:**

Has steroid 11-beta-hydroxylase activity. In addition to this activity, the 18 or 19-hydroxylation of steroids and the aromatization of androstendione to estrone have also been ascribed to cytochrome P450 XIB.

**Gene ID:**

CYP11B1/CYP11B2

**Uniprot**

P15538/P19099

**Synonyms:**

C11B1; CPN1; CYP11B; CYP11B1; CYPXIB1

**Immunogen:**


Synthesized peptide derived from C-terminal of human Cytochrome P450 11B1/2.

**Storage:**

Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---



Western blot analysis of extracts from COLO205 cells, using Cytochrome P450 11B1/2 antibody.