

HSPBAP1 Antibody



PACO22746

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

May play a role in cellular stress response By similarity.

Gene ID:

HSPBAP1

Uniprot

Q96EW2

Synonyms:

HSPB1-associated protein 1; 27 KdA heat shock protein-associated protein 1; Protein associated with small stress protein 1; HSPBAP1; PASS1


Immunogen:

Synthesized peptide derived from internal of human HBAP1.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from HeLa cells and 293 cells, using HBAP1 antibody.