
Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Catalyzes the isomerization of prostaglandin H2 to prostacyclin (= prostaglandin I2).

Gene ID:

PTGIS

Uniprot

Q16647

Synonyms:

CYP8; CYP8A1; cytochrome P450; family 8; subfamily A


Immunogen:

Synthesized peptide derived from Internal of human PTGIS.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from K562 cells, using PTGIS antibody.