

# CYP4A11/CYP4A22 Antibody



PACO22779

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:3000

**Protein Background:**

Catalyzes the omega- and (omega-1)-hydroxylation of various fatty acids such as laurate, myristate and palmitate. Has little activity toward prostaglandins A1 and E1. Oxidizes arachidonic acid to 20-hydroxyeicosatetraenoic acid (20-HETE).

**Gene ID:**

CYP4A11/CYP4A22

**Uniprot**

Q02928/Q5TCH4

**Synonyms:**

CYP4A11 protein;

**Immunogen:**


Synthesized peptide derived from internal of human CYTOCHROME P450 4A11/22.

**Storage:**

Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---



Western blot analysis of extracts from COLO205, HuvEc and HeLa cells, using CYTOCHROME P450 4A11/22 antibody.