

EPHB1 Antibody



PACO22797

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:**Gene ID:**

EPHB1

Uniprot

P54762

Synonyms:

ELK; EPB1; Eph receptor B1; Ephrin type-B receptor 1 precursor; EPHT2


Immunogen:

Synthesized peptide derived from C-terminal of human EPHB1.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from COS7 cells, using EPHB1 antibody.