
Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Prokineticins are secreted proteins that can promote angiogenesis and induce strong gastrointestinal smooth muscle contraction. The protein encoded by this gene is an integral membrane protein and G protein-coupled receptor for prokineticins. The encoded protein is similar in sequence to GPR73, another G protein-coupled receptor for prokineticins.

Gene ID:

PROKR2

Uniprot

Q8NFJ6

Synonyms:

HH3; KAL3; PKR2; GPRg2; GPR73b; GPR73L1; dJ680N4.3


Immunogen:

Synthetic peptide corresponding to residues near the N terminal of human prokineticin receptor 2.

Storage:

Rabbit IgG in pH7.3 PBS, 0.05% NaN₃, 50% Glycerol.

Product Images



Gel: 10%SDS-PAGE Lysate: 40ug 823 cell Primary antibody: 1/200 dilution Secondary antibody dilution: 1/8000 Exposure time: 20 minutes.