

JUN (Ab-239) Antibody



PACO22869

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000,
IHC:1:50-1:200

Protein Background:

Transcription factor that recognizes and binds to the enhancer heptamer motif 5'-TGA[CG]TCA-3'.

Gene ID:

JUN

Uniprot

P05412

Synonyms:

AH119; AP1; Jun A; c-Jun; p39


Immunogen:

Peptide sequence around aa.237~241 (G-E-T-P-P) derived from Human c-Jun.


Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from 293, HT29 and HepG2 cells using c-Jun(Ab-239) Antibody.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using c-Jun(Ab-239) Antibody(left) or the same antibody preincubated with blocking peptide(right).