

# ZAP70 (Ab-319) Antibody



PACO22935

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:1000,  
IHC:1:50-1:200

**Protein Background:****Gene ID:**

ZAP70

**Uniprot**

P43403

**Synonyms:**

70 kDa zeta-associated protein; SRK; Syk-related tyrosine kinase; ZA70; ZAP-70

**Immunogen:**


Peptide sequence around aa.317~321 (S-P-Y-S-D) derived from Human ZAP70.

**Storage:**


Supplied at 1.0mg/mL in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---



Western blot analysis of extracts from HL60 cells using Zap-70(Ab-319) Antibody.



Immunohistochemical analysis of paraffin-embedded human tonsil tissue using Zap-70(Ab-319) Antibody(left) or the same antibody preincubated with blocking peptide(right).