

VASP (Ab-157) Antibody



PACO22951

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000,
IF:1:100-1:200

Protein Background:**Gene ID:**

VASP

Uniprot

P50552


Synonyms:**Immunogen:**

Peptide sequence around aa.155~159 (R-V-S-N-A) derived from Human VASP.


Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4,
150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from 293 and C6 cells using VASP(Ab-157) Antibody.



Immunofluorescence staining of methanol-fixed HeLa cells using VASP(Ab-157) Antibody.