

# HSPB1 (Ab-78) Antibody



PACO22957

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:1000,  
IHC:1:50-1:200

**Protein Background:****Gene ID:**

HSPB1

**Uniprot**

P04792

**Synonyms:**

Heat shock 27 kDa protein; HSP 27; Stress-responsive protein 27

**Immunogen:**


Peptide sequence around aa.76~80 (A-L-S-R-Q) derived from Human HSP27.

**Storage:**


Supplied at 1.0mg/mL in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---



Western blot analysis of extracts from HeLa cells using HSP27(Ab-78) Antibody.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using HSP27(Ab-78) Antibody(left) or the same antibody preincubated with blocking peptide(right).