

ARNTL Antibody



PACO22992

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:

Gene ID:

ARNTL

Uniprot

O00327

Synonyms:

MOP3; BMAL1c; ARNTL; TIC;


Immunogen:

Peptide sequence around aa.569~573 (S-S-S-I-L) derived from Human BMAL1.


Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from mouse brain tissue using BMAL1 antibody.



Western blot analysis of extracts from HL60, THP-1.