

POU2F1 Antibody



PACO23013

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:1000,
IF:1:100-1:200

Protein Background:**Gene ID:**

POU2F1

Uniprot

P14859

Synonyms:

1-Oct; OTF1; oct-1B


Immunogen:

Peptide sequence around aa.270~274(T-P-K-R-I) derived from Human Oct-1.


Storage:

Supplied at 1.0mg/mL in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from JK, HepG2 and 293 cells using Oct-1 Antibody.



Immunofluorescence staining of methanol-fixed NCI-h446 cells showing cytoplasmic staining using Oct-1 Antibody.