

PACO23031

Product Information

Size:

100ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:50-1:100,
IF:1:100-1:500

Protein Background:**Gene ID:**

ALK

Uniprot

Q9UM73

Synonyms:

ALK tyrosine kinase receptor precursor; Anaplastic lymphoma kinase; CD246; EC 2.7.10.1; kinase ALK

Immunogen:


Synthesized peptide derived from Internal of human ALK.

Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images

Immunohistochemistry analysis of paraffin-embedded human brain tissue using ALK antibody.



Immunofluorescence analysis of HeLa cells, using ALK antibody.

