

# DAPK3 (Ab-265) Antibody



PACO23033

---

## Product Information

**Size:**

100ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:3000,  
IHC:1:50-1:100

**Protein Background:****Gene ID:**

DAPK3

**Uniprot**

O43293

**Synonyms:**

DAP kinase 3; DAP- like kinase; Death-associated protein kinase 3; Dlk; EC 2.7.11.1

**Immunogen:**


Synthesized non-phosphopeptide derived from human DAPK3 around the phosphorylation site of threonine 265 (R-M-T(p)-I-A).

**Storage:**

Rabbit IgG in phosphate buffered saline (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Product Images

---



Western blot analysis of extracts from HuvEc cells, using DAPK3 (Ab-265) antibody.

Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue using DAPK3 (Ab-265) antibody.

