

PTK2B (Ab-579) Antibody



PACO23073

Product Information

Size:

100ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:**Gene ID:**

PTK2B

Uniprot

Q14289

Synonyms:

CADTK; CAK beta; Calcium-dependent tyrosine kinase; Cell adhesion kinase beta; EC 2.7.10.2


Immunogen:

Synthesized non-phosphopeptide derived from human PYK2 around the phosphorylation site of tyrosine 579 (E-D-Y(p)-Y-K).


Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from 3T3 cells, using PYK2 (Ab-579) antibody.



Western blot analysis of extracts from 3T3 cells (Lane 2), using PYK2 (Ab-579) antibody. The lane on the left is treated with synthesized peptide.