

KRT19 Antibody



PACO23119

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000,
IHC:1:50-1:100

Protein Background:

Involved in the organization of myofibers. Together with KRT8, helps to link the contractile apparatus to dystrophin at the costameres of striated muscle.

Gene ID:

KRT19

Uniprot

P08727

Synonyms:

Keratin type I cytoskeletal 19; Cytokeratin-19; CK-19; K19;


Immunogen:

Synthesized peptide derived from human Keratin 19.


Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from LOVO cells, using Keratin 19 antibody.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using Keratin 19 antibody.