

TGFB3 Antibody



PACO23138

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000,
IHC:1:50-1:100

Protein Background:

Involved in embryogenesis and cell differentiation.

Gene ID:

TGFB3

Uniprot

P10600

Synonyms:

TRANSFORMING GROWTH FACTOR BETA-3;


Immunogen:

Synthesized peptide derived from human TGF beta 3.


Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from COS7 cells, using TGF beta 3 antibody.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using TGF beta 3 antibody.