

PIAS2 Antibody



PACO23145

Product Information

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:3000

Protein Background:**Gene ID:**

PIAS2

Uniprot

O75928

Synonyms:

ARIP3; Androgen receptor-interacting protein 3; DAB2-interacting protein; DIP; Miz1


Immunogen:

Synthesized peptide derived from Internal of human PIAS2.


Storage:

Rabbit IgG in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Product Images



Western blot analysis of extracts from COS7 cells, using PIAS2 antibody.



Western blot analysis of extracts from HuvEc cells (Lane 2), JK cells (Lane 3) and cos-7 cells (Lane 4), using PIAS2 antibody. The lane on the left is treated with synthesized peptide.